

Allergy File Specification

Field Name	Field Description	Example Output
practice_id	Practice ID	1111
create_timestamp	Blank	
modify_timestamp	Date allergy was last modified	2023-07-25 12:04:00
allergy_date	Date of allergy	2023-07-25
allergy_uuid	UUID of the allergy	ba50c68e-0207-49d1-b3e0-d59bb13cae0b
patient_mrn	Patient MRN	MRN.1111.fc4f8851-3762-4e82-99e0-7a0e1551a202
MEDENT_id	Patient's MEDENT account number	6259
type	Type field from Allergy Detail screen	Adverse Reaction
description	Allergy field from Allergy Detail screen (This is the name of the allergy)	Amoxicillin
reaction_code	SNOMED code for the reaction listed, if any. If multiple exist, they are comma separated.	81555016, 2576576015
reaction_desc	Reaction description from Allergy Detail screen. If multiple exist, they are comma separated.	Hives, Itching
reaction_date	Reaction date in the allergy detail	2014-08-12
snomed_code	SNOMED code of the allergy, if one exists	
unii_code	UNII code of the allergy, if one exists	
rxnorm_code	RxNorm associated with Allergy, if one exists	723
onset_date	Active Date from the Allergy Detail screen	2014-08-12
inactive_date	Inactive Date from the Allergy Detail screen	
severity	Severity from the Allergy Detail screen	Moderate
criticality	Criticality from the allergy detail	Unable to Assess Risk
verification_status	Verification Status from allergy detail (Confirmed, Unconfirmed, Refuted)	Unconfirmed
allergy_comments	Comments field within the Allergy detail	

Notes: If a patient has no known drug allergies (NKDA) indicated in their chart this will be included in the report sending NKDA in both the type and description fields, and the date the NKDA entry was added will send in the onset_date field.