

## CareTeam and Resources File Specification

Field Name	Field Description	Example Output
practice_id	Practice ID	1111
patient_mrn	Patient MRN	MRN.1111.fc4f8851-3762-4e82-99e0-7a0e1551a202
MEDENT_id	Patient's MEDENT account number	6259
provider_id	Doctor # in MEDENT for the entry with an "R" prefix if pulling the referring doctor master file number	3
name	Provider name	John Smith, MD
care_team	Care Team Member	Yes
community_resource	Community Resource	No
notes	Include all notes in this area	