

Document Listing File Specification

Field Name	Field Definition	Example Output
practice_id	Practice ID	1111
create_timestamp	Date/time the document was created	2023-07-25 12:04:00
modify_timestamp	Last modified date of the document	2023-07-25 12:04:00
document_uuid	UUID of the document	a943553a-7b6c-49b9-8155-4659f635f67a
patient_mrn	Patient MRN	MRN.1111.fc4f8851-3762-4e82-99e0-7a0e1551a202
MEDENT_id	Patient's MEDENT account number	6259
encounter_id	e-Superbill# associated with the same date/dr as document	6960
provider_id	Doctor # in MEDENT from document route screen	4
document_code	Document code from patient chart	PCS-ACT
document_name	Sends the document name from patient chart	Asthma Control Test
document_date	Date of the document from route box	2023-07-25
loinc_code	LOINC code attached to the document in the master file setup	82674-3

All document codes with the exception of documents assigned to chart button "NOTES" or any of the following codes: OBGNU-NURS, UROLNU-NOT, NOB, NOF, NFP, NFF, NU1, NU2, PT NOTE*, OT NOTE*, SEND, CCD-INCOMI, CCD-OUTGOI, CCD-PORTAL, ORDAIN, ORDAOUT, UD-IN, UD-INCOMIN, UD-OUT, UD-OUTGOIN, and any document code beginning with NOTE*