

Encounter File Specification

Field Name	Field Description	Example Output
practice_id	Practice ID	1111
create_timestamp	Date/time the encounter was created	2023-03-21 11:26:00
modify_timestamp	Last modified date of the e-Superbill	2023-10-04 10:03:00
encounter_id	E-Superbill number	2339
patient_mrn	Patient MRN	MRN.1111.fc4f8851-3762-4e82-99e0-7a0e1551a202
MEDENT_id	Patient's MEDENT account number	6259
provider_id	Provider number on E-Superbill	1224
provider_npi	Provider NPI from Doctor Master file	1926412317
supervising_id	Supervising Provider number on E-Superbill	
supervising_npi	Supervising Provider NPI from Doctor Master file	
location_id	Location number on the E-Superbill	2
encounter_date	Date of encounter from E-Superbill	2023-03-21
encounter_time	Time of encounter from E-Superbill	11:26:00
encounter_reason	Free Text note in E-Superbill	

Note: E-Superbills with status of Closed or PatSeen are included, as well as any selected in "Additional Ebill Status".