

Exchanged Patient Level Files

Field Name	Field Definition	Example Output
practice_id	Practice ID	1111
create_timestamp	Date/time the document was created	2023-03-21 11:26:00
modify_timestamp	Last modified date of the document	2023-03-21 11:26:00
document_uuid	UUID of the document	4e7f767d-11de-4756-8929-fca0a01772de
patient_mrn	Patient MRN	MRN.1111.fc4f8851-3762-4e82-99e0-7a0e1551a202
MEDENT_id	Patient's MEDENT account number	6259
document_code	Document code from patient chart	SEND
document_name	Document name from patient chart	Sent to ACO
date_of_disclosure	Date of the document from route box	2023-03-21
sent_to	Sent to field within a Send Package	ACO
purpose	Purpose for which the Send package was created	ACO
content_exported	Date and Description of data included in send package. Multiple files exported are listed with the date then filename each separated by semi colons	06/04/23 .Note, Burn; 05/19/23 .Note, Dermatitis Dyshidrotic; 02/08/23 .Note, Burn
loinc_code	LOINC code attached to the document in the master file setup	

This file should include any of the following document codes as well as items linked to receiving data in the Structured Document Setup area: SEND, CCD-INCOMI, CCD-OUTGOI, CCD-PORTAL, CCD-API, QRDAIN, QRDAOUT, UD-IN, UD-INCOMIN, UD-OUT, UD-OUTGOIN