

MedicalHistory File Specification

Field Name	Field Definition	Example Output
practice_id	Practice ID	1111
create_timestamp	Create date and time of the entry	2023-02-07 12:04:00
modify_timestamp	Last modified date and time	2023-11-28 08:55:00
medical_history_uuid	UUID of the entry	5661f374-2193-4eb8-abaf-2919c8f346d1
patient_mrn	Patient MRN	MRN.1111.b2aafb65-0f66-4dad-9dc9-344a28f3cd36
MEDENT_id	Patient's MEDENT account number	1251
description	Description of the medical history item	Abdominal Aortic Aneurysm
history_date	Date of the medical history item	2023-02-07
snomed_ct_code	SNOMED code of the medical history item	87522002
medical_history_text	comments section	