

## Medication File Specification

Field Name	Field Definition	Example Output
practice_id	Practice ID	1111
create_timestamp	Create date and time of the entry	2021-06-29 00:00:00
modify_timestamp	Last modified date and time	2022-12-02 09:55:00
medication_uuid	UUID of the entry	cabd15f4-01d3-4ba3-b085-2a2de607e7e8
patient_mrn	Patient MRN	MRN.1111.82440cab-3f8f-4d7b-9e6f-96512016479d
MEDENT_id	Patient's MEDENT account number	35
encounter_id	e-SuperBill# for same date/dr/loc as Medication entry	1259
provider_id	Provider ID in MEDENT	1
other_prov	1 if the entry has been prescribed by another provider, 0 otherwise	0
loc_id	Location number	2
ndc_code	ndc code	591040801
rxnorm_code	rxnorm code	314077
dose	dose (e.g. 17.5-3.13-1.6)	20
units	units (e.g. GM/177ML)	mg
type	type (e.g. Oral Solution)	Tablets
quantity	quantity	60
route	route	Oral
rx_method	how the prescription was sent (eRx, Fax, Print)	eRx
start_date	prescribed date	2021-04-29
stop_date	D/C date if one exists	2022-12-02
disc_reason	D/C reason, if one exists	
disc_comment	Comments from the D/C area of a medication	
dispensed_date	Dispensed date	2021-04-29
dea_schedule	Dea schedule	
rx_refills	# of refills noted in the med	3
date_reordered	Date med was reordered	
form	Form from Medication name in prescription	Tablets
sig_description	Description of the SIG for medication	1 tab by mouth twice a day as directed
mdd	MDD field in the medication detail for maximum daily dosage	1
supply_days	Days Supply field in the prescription record	
signed_by	Provider name that signed the prescription	
signed_timestamp	Date the provider signed the prescription	
admin_route	Route from Administered Values, when administered in office	Oral
admin_site	Site from Administered Values, when administered in office	
admin_start	Start date from Administered Values, when administered in office	04/29/21
admin_stop	Stop date from Administered Values, when administered in office	
admin_lot	Lot# from Administered Values, when administered in office	

admin_expdate	Expired Date from Administered Values, when administered in office	
medication_name	Description of the medication, as shown in MEDENT	Lisinopril
pharmacy_sent	Pharmacy name the prescription was sent to	Walgreens
dx	Diagnosis attached to the medication. If multiple, they should be comma delimited.	
as_needed	Yes or No if the As Needed box is checked in the medication detail.	No
sample_info	Lot# and Expiration date, if marked as a sample.	
pt_education	Yes or No if the box for include patient education is checked	Yes