

OBEpisode File Specification

| Field Name | Field Definition |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| practice_id | Practice ID |
| create_timestamp | Create date of the Antepartum Record document |
| modify_timestamp | Last modified date of the document |
| ob_episode_uuid | UUID of Antepartum Record document |
| patient_mrn | Patient MRN |
| MEDENT_id | Patient's MEDENT account number |
| ob_provider_id | Provider# from route box of Antepartum Record document |
| pregnancy_number | Gravidity entry from PFSH document |
| initial_exam_date | Date from OBG-ANTEPARTUM FIRST PRENATAL VISIT date component. Complete search of data is as follows: <ol style="list-style-type: none"> 1. Search date in OBG-ANTEPARTUM FIRST PRENATAL VISIT 2. If nothing, search OBG-INITAL EXAM PREFORMED 3. If nothing, search OBG-INITIAL EXAM PERFORMED 4. If nothing, search Initial Exam Date from the EDD Calculator |
| initial_gest_age | Content of component OBG-GESTATIONAL AGE WEEKS ONLY if the OBG-ANTEPARTUM FIRST PRENATAL VISIT date has been entered. If nothing is found, the search continues as follows: <ol style="list-style-type: none"> 1. OBG-GESTATIONAL AGE WEEKS ONLY and pull ONLY if the OBG-ANTEPARTUM FIRST PRENATAL VISIT includes a date 2. Next, checks GA field next to the Initial Exam Date field in the EDD Calculator |
| initial_exam_with | OBG-ANTEPF document is searched for the following data: <ol style="list-style-type: none"> 1. OBG-ANTEPARTUM FIRST PRENATAL VISIT component, if there is something there, this field populates 'OTHER' 2. OBG-INITAL EXAM PREFORMED component, if there is a date, populate with 'CENTER' 3. Initial Exam Date field in the EDD Calculator, if there is a date, populate with 'CENTER' |
| edd | Expected Delivery Date is found in the following order until a date is identified: <ol style="list-style-type: none"> 1. Final EDD from the EDD Calculator 2. OB-EDD FINAL EDD(date component) in the OBG-ANTEPF document 3. EDD by 2nd Ultrasound from EDD calculator 4. EDD by 1st Ultrasound from EDD Calculator 5. EDD by Initial Exam from EDD Calculator 6. EDD by Conception from EDD Calculator 7. EDD by LMP from EDD Calculator |
| ob_risk | OBG-ANTEPF document is searched to capture selection component OBSTETRIC PROBLEM LIST; one entry will be created for each item noted in chart/document component. |

Includes the following Antepartum Record documents: OBG-ANTEPF; OBG-ANTEPM; OBG-ANTEPT