

Order File Specification

| Field Name | Field Definition | Example Output |
|-------------------------|--|---|
| practice_id | Practice ID | 1111 |
| create_timestamp | Create date and time of the entry | 2017-03-23 16:26:00 |
| modify_timestamp | Last modified date and time | 2023-07-27 14:16:00 |
| order_uuid | UUID of the entry | b793e8a9-7dbf-46ad-8bb6-9f40f09e45f0 |
| encounter_id | e-Superbill# if found on same date/dr/loc as the order | 259 |
| patient_mrn | Patient MRN | MRN.1111.1df9ea01-5c42-4414-8e2f-556fcd4c07c6 |
| MEDENT_id | Patient's MEDENT account number | 51 |
| location_id | Location in the DI/Order/Lab | 1 |
| provider_id | Provider ID in MEDENT for the PDr in the order | 1 |
| order_provider_id | Provider ID in MEDENT for the Odr in the order | 1 |
| supervising_provider_id | Provider ID in MEDENT for the SDr in the order | 245 |
| referring_provider_id | Provider ID in MEDENT for the RDr in the order | R35 |
| facility | Facility name for the DI/Order or lab company name | Patients Choice |
| fasting | Contents of the Fasting field within a lab order | |
| type | Diagnostic Image/General Order/Lab Order based on where the data is coming from in the chart | Diagnostic Image |
| loinc_code | Loinc code attached to the order | 38268-9 |
| order_code | House code ordered | 77081 |
| order_desc | Order description | Dexa Bone Density Study Appendicular Skeleton |
| ordered_date | Ordered date | 2017-03-23 |
| scheduled_date | Scheduled date | |
| acquired_date | Acquired date | |

| | | |
|---------------|--|----------|
| received_date | Received complete/received part date | |
| closed_date | Closed date | |
| declined_date | Declined date | |
| icd10_codes | ICD10 codes attached to the order (can be up to 3 included). Comma separated if multiple | M80.00xD |
| cpt_codes | CPT codes attached to the order (can be up to 3 included). Comma separated if multiple | 77081 |
| priority | Priority rating in detail screen | Routine |
| auc_status | AUC status field | Adheres |
| cdsm_tool | CDSM tool indicated | NDSC |
| auc_code | AUC Gcode indicated | G1004 |
| auc_modifier | AUC mofifier | ME |
| accession_id | Accession/Lab# field from lab detail screen, blank if not documented | 9-9 |