

## PatientPlan File Specification

Field Name	Field Description	Example Output
practice_id	Practice ID	1111
create_timestamp	Date/time the insurance was added to the patient	2016-05-13 09:17:00
patient_mrn	Patient MRN	MRN.1111.1df9ea01-5c42-4414-8e2f-556fcd4c07c6
MEDENT_id	Patient's MEDENT account number	51
plan_id	Insurance ID# from patient insurance screen. This is a numerical identifier attached to the insurance in the master files that is different for each MEDENT practice.	219
policy_number	Subscriber ID# field in patient's insurance member detail screen.	ZXY123456789
policy_start_date	Effective Date from patient's insurance member detail screen	2016-05-13
policy_end_date	Expired Date from patient's insurance member detail screen	
rank	Billing Order from patient's insurance member detail screen	1