

ProblemList File Specification

Field Name	Field Definition	Example Output
practice_id	Practice ID	1111
create_timestamp	Create date and time of problem entry	2021-06-02 09:08:00
modify_timestamp	Last modified date and time	2021-06-02 09:08:00
problem_uuid	UUID of the entry	54372244-68c5-4c86-9d78-efe28679d4cc
encounter_id	e-Superbill# for same date/dr/loc of the problem list entry	2546
patient_mrn	Patient MRN	MRN.1111.1df9ea01-5c42-4414-8e2f-556fcd4c07c6
MEDENT_id	Patient's MEDENT account number	51
provider_id	Provider's MEDENT ID	1
other_prov	1 if the other provider is checked in the entry, 0 otherwise	0
snomed_ct_code	SNOMED code associated with the problem list entry	82959004
description	Description of the SNOMED code associated with the problem list entry	Dementia paralytica juvenilis
onset_date	Onset date from problem list entry	2021-01-15
resolved_date	Status date field of problem list entry when status is resolved	
chronic_ind	1 if Chronic is checked, 0 otherwise	0
rank	Rank of the diagnosis in the problem list	18
reviewed_date	Last reviewed date	