

Progress Note Listing File Specification

Field Name	Field Definition	Example Output
practice_id	Practice ID	1111
create_timestamp	Date/time the document was created	2021-06-29 00:00:00
modify_timestamp	Last modified date of the progress note	
document_uuid	UUID of the document	acdbfbff-8841-4790-ae11-0224d00d2dbd
patient_mrn	Patient MRN	MRN.1111.1df9ea01-5c42-4414-8e2f-556fcd4c07c6
MEDENT_id	Patient's MEDENT account number	51
encounter_id	e-Superbill# when on same date/dr as progress note	1313
provider_id	Dr # in MEDENT from document route screen	1
super_dr	Supervising Dr # from document route screen	
loc_id	Locations ID from the document route screen	
dx	Diagnosis code(s) in the document route screen	Z71.3, E08.3212, E13.21
document_code	Document code from patient chart	NOTE
document_name	Document name from patient chart	.Note, .Exam_Initial
document_date	Date of the document from route box	2021-06-27
loinc_code	LOINC code attached to the document in the master file setup	11506-3
doctor_signature	Yes/No if there is a Doctor Signature within the document	Yes

This file includes all documents assigned to the chart button "NOTES" or with any of the following codes: OBGNU-NURS, UROLNU-NOT, NOB, NOF, NFP, NFF, NU1, NU2, PT NOTE*, OT NOTE* and any document code beginning with NOTE*