

Referral File Specification

Field Name	Field Definition	Example Output
practice_id	Practice ID	1111
create_timestamp	Create date time	2021-12-14 00:00:00
modify_timestamp	Last modified date time	
referral_id	Unique ID of the referral	141
patient_mrn	Patient MRN	MRN.1111.1df9ea01-5c42-4414-8e2f-556fcd4c07c6
MEDENT_id	Patient's MEDENT account number	51
encounter_id	e-Superbill# for the encounter on same date/provider/location as the referral	1537
referral_number	Referral# in the insurance section	
from_location_id	Location from within the referral	1
from_provider_id	Dr # in the RDr field	R1
to_provider_id	Dr# in the Ref to Dr field	R30
to_provider	Dr name in the Ref to Dr Field	Jones, Jessica MD
ref_type_code	Referral Type SNOMED code	453641000124107
ref_type_desc	Referral Type description	Referral to specialist pharmacist
ref_reason	Reason for Referral from Referral section	
ordered_date	created date	2021-12-14
scheduled_date	scheduled date	
sent_date	Sent date	
appointment_date	date in the appointment section	
received_date	received complete date	
closed_date	closed date	
declined_date	declined date	
referral_link	Sends Y if anything is listed in Links sections of referral detail screen, N if not	N
visits	Visits	Unlimited
effective_date	effective date	2021-12-14
expired_date	expired date	2022-12-14