

## Result File Specification

Field Name	Field Definition	Example Output
practice_id	Sends practice's ID	1111
create_timestamp	create date time	2023-07-27 14:16:00
modify_timestamp	Last modified date time	2023-07-27 14:16:00
result_id	UUID for the resulted test. This could be the same as the order UUID, but often is not. (The UUID on the order is likely from a Lab Panel and UUID for the resulted test is on the individual lab test within the panel.)	f9a38fde-6eda-40ff-be93-eaf2d4921ddc
order_uuid	UUID for the ordered test	f9a38fde-6eda-40ff-be93-eaf2d4921ddc
patient_mrn	Patient MRN	MRN.1111.1df9ea01-5c42-4414-8e2f-556fcd4c07c6
MEDENT_id	Patient's MEDENT account number	51
encounter_id	e-Superbill# if found on same date/dr/loc as the order	
type	Diagnostic Image/General Order/Lab Order based on the area of the chart the data came from	Lab Order
facility	Facility name for the DI/Order or lab company name	Patients Choice
location_id	Location in the DI/Order/Lab	1
loinc_code	Loinc code attached	4548-4
order_code	House code ordered	A1C W/RFX
description	Lab description	Hemoglobin A1c W/RFX Glycomark
result	Result field	10.1
high_low	H/L	
reference_range	Reference range	
priority	Priority	Routine
units	Units field	
status	Last status	RcvdComp
notes	Result notes	
ordered_date	Ordered date	2017-01-01
scheduled_date	Scheduled date	
acquired_date	Acquired date	2017-01-05
received_date	Received complete date	2017-01-05
reviewed_date	Reviewed complete date	
closed_date	Closed date	
declined_date	Declined date	
provider_id	Provider ID in MEDENT for the PDr in the order	1
order_provider_id	Provider ID in MEDENT for the Odr in the order	1
supervising_provider_id	Provider ID in MEDENT for the SDr in the order	1
referring_provider_id	Provider ID in MEDENT for the RDr in the order	
icd10_codes	ICD10 codes attached (can be up to 3 included). Comma separated if multiple	E11.9
cpt_codes	CPT codes attached (can be up to 3 included). Comma separated if multiple	

auc_status	AUC status field	
cdsm_tool	CDSM tool indicated	
auc_code	AUC Gcode indicated	
auc_modifier	AUC modifier	
accession_id	Sends Accession/Lab# field from lab detail screen, blank if not documented	