

## SocialHistory File Specification

Field Name	Field Definition	Example Output
practice_id	Practice ID	1111
create_timestamp	Create date time for the entry	2023-11-29 07:22:00
modify_timestamp	Last modified date time for the entry	2023-11-29 07:22:00
social_history_uuid	UUID of the Social History entry	71835045-7d2d-4933-8fb6-5750b2de4588
patient_mrn	Patient MRN	MRN.1111.1df9ea01-5c42-4414-8e2f-556fcd4c07c6
MEDENT_id	Patient MEDENT account number	51
type	Description of the social history item	Lives With: Boyfriend
history_date	Date associated with the social history item	2023-11-29
snomed_ct_code	SNOMED code associated with social history item	307104007
social_history_text	comments section	