

Vital File Specification

Field Name	Field Definition	Example Output
practice_id	Practice ID	1111
create_timestamp	create date time	12/14/2021 14:32
modify_timestamp	Last modified date time	
vitals_uuid	UUID of the individual vital entry	f125c337-00d1-4722-a16c-526d4965116e
patient_mrn	Sends MRN without "MRN." prefix	MRN.1111.1df9ea01-5c42-4414-8e2f-556fcd4c07c6
MEDENT_id	Sends patient's MEDENT account number	51
encounter_id	Sends E-Superbill # with status "PatSeen" or "Closed" when same date and doctor as vital entry	1537
provider_id	provider # in MEDENT (R before the provider # if a provider from the referring master file has been selected)	1
other_provider	1 if other provider checked, 0 otherwise	0
remote_monitored_result	1 if result is noted as being remote monitored and 0 otherwise.	0
vital_name	Description of the vital (e.g. Blood Pressure Recheck)	Blood Pressure
vital_result	Result of the vital	130/85
vital_unit	Units (e.g. mmHg)	mmHg
vital_comment	Comments noted, if any	
vital_loinc	Loinc for the vital entry, if there is one	55284-4