

(g)(10) Standardized API - MEDENT VERSION 23.5 R4**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

CarePlan

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the patient CarePlan file record.
text	Narrative	FALSE	Name of the Care Plan and concern being addressed.
status	string	FALSE	Status of the Care Plan. Values include active, completed, or errored. o Entered-in-error is supported, but errored records will be filtered from responses.
intent	string	FALSE	Value will be plan.
category*	CodeableConcept	FALSE	Value will be assess-plan.
title	string	TRUE	Name of the Care Plan.
subject	Reference	FALSE	References the Patient Resource.
activity*	Activity	FALSE	What is to occur as part of the care plan.

Activity

Name	Type	Is Optional?	Description
detail	Detail	TRUE	Details for the care plan.

Detail

Name	Type	Is Optional?	Description
kind	string	TRUE	Will be "Task"
status	string	FALSE	Value will be indicated as Unknown.