

**(g)(10) Standardized API - MEDENT VERSION 23.7**

**Bundle**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

**CareTeam**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
identifier*	Identifier	Tru	UUID of the Care Team record
status	string	TRUE	Response will include active members
category*	CodeableConcept	TRUE	Value will be "longitudinal".
name	string	TRUE	Value will be "Care Team".
subject	Reference	FALSE	References the Patient resource.
managingOrganization	Reference		References the Organization resource.
participant*	Participant	FALSE	Who is included in the CareTeam

**Participant**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
-------------	-------------	---------------------	--------------------

member	Reference	FALSE	This will reference the Practitioner Resource for Referring or Performing Provider(s), Organization for an Insurance Company or Rolodex entry, or the Patient with RelatedPerson resource
role	CodeableConcept	FALSE	Values may include Healthcare professional (SNOMED 223366009), Services (SNOMED 224930009), Patient (SNOMED 116154003), or Caregiver (SNOMED 133932002)