

(g)(10) Standardized API - MEDENT VERSION 23.7**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Condition

Name	Type	Is Optional?	Description
identifier	Identifier	TRUE	UUID of the Problem List entry
clinicalStatus	CodeableConcept	TRUE	Values include Active, Resolved, or Inactive.
verificationStatus	CodeableConcept	TRUE	Active and Resolved statuses will be sent as Confirmed.
category*	CodeableConcept	FALSE	Value will be "problem-list-item" or will be "health-concern" if entry is linked to the Concern of an Interactive Care Plan.
category:screening-assessment	CodeableConcept	TRUE	For a health concern entry, this is will be shared if the SNOMED is identified as a social determinant of health code "sdoh".
code	CodeableConcept	FALSE	Coding for the problem list entry. This will be SNOMED coded.
subject	Reference	FALSE	References the Patient Resource
onset	Period	TRUE	Onset period of the condition
abatement	Period	TRUE	Abatement of the condition if inactive or resolved
note*	Annotation	TRUE	Note/Comment associated with the record

recordedDate	DateTime	TRUE	If the UUID for the condition has a create date that is after or more recent than than the documented Onset date in the problem list detail screen, uses the date in the detail screen. Otherwise we'll use the Created Date of the UUID. We will also use the Created Date of the UUID if the Onset is set to 00/00/0000.
--------------	----------	------	--