

(g)(10) Standardized API - MEDENT VERSION 23.7

Bundle

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Condition

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the Diagnosis record
clinicalStatus	CodeableConcept	TRUE	Value will be Active.
verificationStatus	CodeableConcept	TRUE	Value will be confirmed.
category*	CodeableConcept	FALSE	Value will be "encounter-diagnosis"
code	CodeableConcept	FALSE	Coding for the problem list entry. This will be ICD10.
subject	Reference	FALSE	References the Patient resource
encounter	Reference	TRUE	References the Encounter where this diagnosis code was selected.
onset	DateTime	TRUE	Date associated with the diagnosis record
abatement	DateTime	TRUE	Date associated with the diagnosis record
note*	Annotation	TRUE	Note/Comment associated with the diagnostic record