

(g)(10) Standardized API - MEDENT VERSION 23.7**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Coverage

Name	Type	Is Optional?	Description
identifier:memberid	Identifier	TRUE	memberID
status	code	FALSE	This will send "active" if the Insurance Effective Period has not elapsed
type	CodeableConcept	TRUE	Source of Payment Typology from insurance master file entry
subscriberId	string	TRUE	memberID
beneficiary	Reference	FALSE	Reference Patient
relationship	CodeableConcept	FALSE	Subscriber relationship from Insurance Detail
period	Period	TRUE	Insurance Effective Start and Stop Dates
payor	Reference	FALSE	PayerID Organization Resource
class:group	BackboneElement	TRUE	Member Group Information
class:plan	BackboneElement	TRUE	Member Group Information