

**(g)(10) Standardized API - MEDENT VERSION 23.7****Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

**DocumentReference**

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the Document
masterIdentifier	Identifier		Versioned identifier of the Document record.
status	code	FALSE	Status will be "current". MEDENT will share the most current version fo the document.
type	CodeableConcept	FALSE	Value will be "problem-list-item".
subject	Reference	FALSE	Coding for the problem list entry. This will be SNOMED coded.
date	instant	TRUE	When the record id was created
custodian	Organization	TRUE	References the Organization resource
author*	Reference	TRUE	References the Practitioner resource
description	string	TRUE	Name of the document
content*	Content	FALSE	Document information referenced

**Content**

Name	Type	Is Optional?	Description
attachment	Attachment	FALSE	

**Attachment**

Name	Type	Is Optional?	Description
------	------	--------------	-------------

contentType	code	FALSE	Mime content type string. Possible mimetypes: application/word application/excel text/xml application/octet-stream text/pdf application/excel application/pdf
data	base64Binary	TRUE	document converted to PDF, base64 encoded
hash	base64Binary	TRUE	SHA-1 hash created from the document.
title	string	TRUE	Name of the document
creation	dateTime	TRUE	UTC date/time when segment was created
url	url	TRUE	Where the data can be found