

**(g)(10) Standardized API - MEDENT VERSION 23.7****Bundle**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

**Encounter**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
identifier*	Identifier	TRUE	UUID of the Superbill for the Encounter
status	code	FALSE	Values include planned, cancelled, entered-in-error, arrived, in-progress, finished, and unknown depending on the status of the Superbil
statusHistory*		TRUE	Past encounter statuses
class	Coding	FALSE	Value will be "ambulatory"
subject	Reference	FALSE	References the Patient resource
participant*	Participant	TRUE	This is based on the Provider associated with the encounter
statusAppointment	Reference	TRUE	References the Appointment resource.
period	Period	TRUE	Start and End time of the Encounter
diagnosis*	Reference	TRUE	References the Procedure resource.

location*	Reference	TRUE	References the Location resource.
serviceProvider	Reference	TRUE	References the Organization resource.
reasonCode	CodeableConcept	TRUE	CPT codes that are linked to Superbills either Directly [CPT/HCPCS screen will have an eBill number], or Indirectly [Doctor, Location, and WCMV match]
reasonReference	Reference	TRUE	References the ICD10 Dx Codes as a Condition Resource Reference

### Participant

Name	Type	Is Optional?	Description
individual	Reference	TRUE	References the Practitioner resource.
type*	code	TRUE	This can indicate if the participant is the Primary Performer (PPRF), Attending Provider (ATND), or the Referring Provider (REF).
period	Period	TRUE	If linked to an Appt, will include the start and end times