

**(g)(10) Standardized API - MEDENT VERSION 23.5 R4**

**Bundle**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

**Encounter**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
identifier*	Identifier	TRUE	UUID of the Superbill for the Encounter
status	code	FALSE	Values include planned, cancelled, entered-in-error, arrived, in-progress, finished, and unknown depending on the status of the Superbil
statusHistory*		TRUE	Past encounter statuses
class	Coding	FALSE	Value will be "ambulatory"
subject	Reference	FALSE	References the Patient resource
participant*	Participant	TRUE	This is based on the Provider associated with the encounter
statusAppointment	Reference	TRUE	References the Appointment resource.
period	Period	TRUE	Start and End time of the Encounter
diagnosis*	Reference	TRUE	References the Procedure resource.
location*	Reference	TRUE	References the Location resource.
serviceProvider	Reference	TRUE	References the Organization resource.

**Participant**

Name	Type	Is Optional?	Description
individual	Reference	TRUE	References the Practitioner resource.
type*	code	TRUE	This can indicate if the participant is the Primary Performer (PPRF), Attending Provider (ATND), or the Referring Provider (REF).