

(g)(10) Standardized API - MEDENT VERSION 23.7**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Goal

Name	Type	Is Optional?	Description
identifier	Identifier*	TRUE	UUID of the Patient Goal
lifecycleStatus	code	FALSE	Values include proposed, planned, accepted, active, on-hold, completed, cancelled, entered-in-error, rejected. Entered-in-error is supported, but errored records will be filtered from responses.
achievementStatus	CodeableConcept	TRUE	Values include in-progress, improving, worsening, no-change, achieved, sustaining, not-achieved, no-progress, not-attainable. Entered-in-error is supported, but errored records will be filtered from responses.
category*	CodeableConcept	TRUE	Values include nursing, dietary, nursing, behavioral, nursing, safety, behavioral, behavioral, behavioral, and nursing.
priority	CodeableConcept	TRUE	Provider indicated priority of the Goal. Values include high, medium, and low.
description	CodeableConcept	FALSE	Goal title
subject	Reference	FALSE	References the Patient resource.
start	date	TRUE	Start Date
target*	Target	TRUE	Goal milestones
statusDate	date	TRUE	When the goal status took effect
statusReason	string	TRUE	Goal reason
note*	Annotation	TRUE	Comment/note associated with th

outcomeCode*	CodeableConcept	TRUE	Outcome information
--------------	-----------------	------	---------------------

Target

Name	Type	Is Optional?	Description
measure	CodeableConcept	TRUE	Values include nursing, dietary, nursing, behavioral, nursing, safety, behavioral, behavioral, behavioral, and nursing.
detailString	string	TRUE	Description of the milestone
dueDate	string	TRUE	Achieve by date for the goal

