

(g)(10) Standardized API - MEDENT VERSION 23.5 R4

Bundle

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Group

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the Group
active	boolean	TRUE	Indicates if active
type	code	FALSE	Value will be "person".
actual	boolean	FALSE	Actual or Descriptive
name	string	TRUE	Name of the Group
managingEntity	Reference	TRUE	References the Organization resource.