

**(g)(10) Standardized API - MEDENT VERSION 23.5 R4****Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

**Immunization**

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the shot record
status	code	FALSE	Will be "completed" if administered and "not-done" if refused
statusReason	CodeableConcept	TRUE	Will display if the immunization was refused.
vaccineCode	CodeableConcept	FALSE	CVX Vaccine Code from the CDC Immunization file.
patient	Reference	FALSE	References the Patient resource.
occurrence	dateTime	FALSE	When the shot was administered or refused
primarySource	boolean	TRUE	Indicates if administered by a provider at the practice or Other Provider
location	Reference	TRUE	References the Location resource
manufacturer	string	TRUE	Manufacturer of the immunization
lotNumber	string	TRUE	Lot Number
expirationDate	date	TRUE	Expiration Date
site	CodeableConcept	TRUE	Site administered
route	CodeableConcept	TRUE	Route administered
doseQuantity	SimpleQuantity	TRUE	Dosage
performer*	Performer	TRUE	Who administered the shot
note*	Annotation	TRUE	
reaction*	Reaction	TRUE	Reaction(s) to the immunization

**Performer**

Name	Type	Is Optional?	Description
------	------	--------------	-------------

actor	Reference	FALSE	References the Practitioner resource. Provider will be the "Administered By" Provider if not marked as Other Provider. If marked as Other Provider, will send the Provider Number if available.
function	CodeableConcept	TRUE	Role of the Provider. Values include Administering Provider (AP) and Ordering Provider (OP).

### Reaction

Name	Type	Is Optional?	Description
detail	Reference	TRUE	References the Observation resource.