

(g)(10) Standardized API - MEDENT VERSION 23.7

Bundle

Name	Type	Is Optional?
id	string	TRUE
meta	instant	TRUE
type	code	FALSE
total	int64	TRUE
link*	url	TRUE

Location

Name	Type	Is Optional?
identifier	Identifier*	TRUE
status	code	TRUE
name	string	FALSE
mode	code	TRUE
telecom*	ContactPoint	TRUE
address	Address	TRUE
managingOrganization	Reference	TRUE
endpoint*	Reference	TRUE

Address

Name	Type	Is Optional?
use	string	TRUE
line	string	TRUE
type	string	TRUE
address-city	string	TRUE
address-state	string	TRUE
address-postalCode	string	TRUE
address-country	string	TRUE
address-period	Period	TRUE

ContactPoint

Name	Type	Is Optional?
system	string	TRUE
value	string	TRUE
use	string	tru

Description
processId-GeneratedDateTime in UTC format
UTC format timestamp of the lastUpdate instant.
searchset
Total number of resources included in the response.
Link for the included data

Description
Location ID
Value will be active or inactive.
Name of the Location
Value will be instance.
Telecom information
Address information
References the Organization Resource.
References the Endpoint Resource

Description
Value will be "work"
Address Line
Value will be "postal"
Name of the City
Name of the State
Postal Code
Country
Period when in use

Description
Value will be phone or fax.
Phone or Fax #
Value will be work.