

**(g)(10) Standardized API - MEDENT VERSION 23.5 R4**

**Bundle**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

**Location**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
identifier	Identifier*	TRUE	Location ID
status	code	TRUE	Value will be active or inactive.
name	string	FALSE	Name of the Location
mode	code	TRUE	Value will be instance.
telecom*	ContactPoint	TRUE	Telecom information
address	Address	TRUE	Address information
managingOrganization	Reference	TRUE	References the Organization Resource.
endpoint*	Reference	TRUE	References the Endpoint Resource

**Address**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
use	string	TRUE	Value will be "work"
line	string	TRUE	Address Line
type	string	TRUE	Value will be "postal"
address-city	string	TRUE	Name of the City
address-state	string	TRUE	Name of the State
address-postalCode	string	TRUE	Postal Code

address-country	string	TRUE	Country
address-period	Period	TRUE	Period when in use

**ContactPoint**

Name	Type	Is Optional?	Description
system	string	TRUE	Value will be phone or fax.
value	string	TRUE	Phone or Fax #
use	string	tru	Value will be work.