

**(g)(10) Standardized API - MEDENT VERSION 23.7**

**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

**Medication**

Name	Type	Is Optional?	Description
identifier	Identifier	TRUE	Document UUID
type	CodeableConcept	TRUE	Value will be "image"
subject	Reference	TRUE	Patient Resource
content*	CodeableConcept	TRUE	

**content**

Name	Type	Is Optional?	Description
attachment	BackboneElement	FALSE	Document UUID
attachment.ContentType	Code	TRUE	Value will be "image"
attachment.data	base64Binary	TRUE	base64 encoded PDF of the attachment
attachment.hash	base64Binary	TRUE	SHA-1 has created from document PDF
attachment.title	string	TRUE	Document Name in medent
attachment.creation	dateTime	TRUE	UTC date/time that this segment was created

This resource will be shared with DiagnosticReports if there are additional attachment links for a Radiology/Diagnostic Imaging order.