

(g)(10) Standardized API - MEDENT VERSION 23.5 R4**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Medication

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	Medispan NDC Code
code	CodeableConcept	FALSE	Medication code will be RxNorm.
manufacturer	CodeableConcept	TRUE	References the Organization resource. Manufacturer information is provided by the Medispan file.
form	CodeableConcept	TRUE	Form of the medication. Manufacturing and form information is provided by the Medispan file.
ingredient*	Ingredient	TRUE	Ingredient from the Medispan file

Ingredient

Name	Type	Is Optional?	Description
item	CodeableConcept	FALSE	Medication ingredient loaded from the Medispan Ingredient file.
isActive	boolean	TRUE	Indicates if the ingredient is active in the Medication