

**(g)(10) Standardized API - MEDENT VERSION 23.7**

**Bundle**

| Name  | Type    | Is Optional? |
|-------|---------|--------------|
| id    | string  | TRUE         |
| meta  | instant | TRUE         |
| type  | code    | FALSE        |
| total | int64   | TRUE         |
| link* | url     | TRUE         |

| Name                    | Type            | Is Optional? |
|-------------------------|-----------------|--------------|
| identifier*             | Identifier      | TRUE         |
| status                  | code            | FALSE        |
| medication              | Reference       | FALSE        |
| subject                 | Reference       | FALSE        |
| authorizingPrescription | Reference       | TRUE         |
| type                    | CodeableConcept | TRUE         |
| performer               | Reference       | FALSE        |
| whenHandedOver          | dateTime        | TRUE         |
| dosageInstruction*      | Dosage          | TRUE         |
| quantity                | Quantity        | TRUE         |

**Dosage**

| Name     | Type    | Is Optional? |
|----------|---------|--------------|
| text     | string  |              |
| isActive | boolean |              |

| Description   |
|---|
| processId-GeneratedDateTime in UTC format           |
| UTC format timestamp of the lastUpdate instant.     |
| searchset   |
| Total number of resources included in the response. |
| Link for the included data                          |

| Description   |
|---|
| UUID of the Medication  |
| Value will indicate if stopped, active, or unknown if no other status can be identified   |
| References the Medication resource.   |
| References the Patient resource.  |
| References the MedicationRequest  |
| Based on status from MedicationDispense History.<br>If Filled - send Completed<br>Partial Filled - Send In-Progress<br>If NotDispensed - send On-Hold |
| References the Organization Resource for the Pharmacy. This is based on mapping the NCPDP (7 digi   |
| Date from the Dispensed History   |
| SIG   |
| Quantity  |

| Description   |
|---|
| This will include the entire MEDENT medication SIG. |
|   |

t National Pharmacy number).