

(g)(10) Standardized API - MEDENT VERSION 23.5 R4**Bundle**

| Name | Type | Is Optional? | Description |
|-------|---------|--------------|---|
| id | string | TRUE | processId-GeneratedDateTime in UTC format |
| meta | instant | TRUE | UTC format timestamp of the lastUpdate instant. |
| type | code | FALSE | searchset |
| total | int64 | TRUE | Total number of resources included in the response. |
| link* | url | TRUE | Link for the included data |

MedicationRequest

| Name | Type | Is Optional? | Description |
|------------------------|-----------------|--------------|--|
| identifier* | Identifier | TRUE | UUID of the Medication |
| status | code | FALSE | Value will indicate if stopped, active, or unknown if no other status can be identified |
| medication | Reference | FALSE | References the Medication resource. |
| subject | Reference | FALSE | References the Patient resource. |
| supportingInformation* | Reference | TRUE | References the Observation resource. Information on the Most Recent Height and Weight observations will be shared. |
| authoredOn | dateTime | TRUE | Date the Medication was prescribed |
| requester | Reference | FALSE | References the Practitioner resource. |
| note* | Annotation | TRUE | |
| dosageInstruction* | Dosage | TRUE | SIG |
| dispenseRequest | DispenseRequest | TRUE | Medication supply information |

Dosage

| Name | Type | Is Optional? | Description |
|------|------|--------------|-------------|
|------|------|--------------|-------------|

| | | | |
|----------|---------|--|---|
| text | string | | This will include the entire MEDENT medication SIG. |
| isActive | boolean | | |

DispenseRequest

| Name | Type | Is Optional? | Description |
|------------------------|----------------|---------------------|---|
| validity | Period | | Prescribed-Discontinued date period |
| numberOfRepeatsAllowed | unsignedInt | | Number of refills indicated for the medication order. |
| quantity | SimpleQuantity | | Quantity prescribed |
| expectedSupplyDuration | Duration | | Expected Discontinue Days |
| performer | Reference | | References the Organization resource. |