

(g)(10) Standardized API - MEDENT VERSION 23.5 R4

Bundle

Name	Type	Is Optional?
id	string	TRUE
meta	instant	TRUE
type	code	FALSE
total	int64	TRUE
link*	url	TRUE

MedicationRequest

Name	Type	Is Optional?
identifier*	Identifier	TRUE
status	code	FALSE
medication	Reference	FALSE
subject	Reference	FALSE
supportingInformation*	Reference	TRUE
authoredOn	dateTime	TRUE
requester	Reference	FALSE
note*	Annotation	TRUE
dosageInstruction*	Dosage	TRUE
dispenseRequest	DispenseRequest	TRUE

Dosage

Name	Type	Is Optional?
text	string	
isActive	boolean	

DispenseRequest

Name	Type	Is Optional?
validity	Period	
numberOfRepeatsAllowed	unsignedInt	
quantity	SimpleQuantity	
expectedSupplyDuration	Duration	
performer	Reference	

Description
processId-GeneratedDateTime in UTC format
UTC format timestamp of the lastUpdate instant.
searchset
Total number of resources included in the response.
Link for the included data

Description
UUID of the Medication
Value will indicate if stopped, active, or unknown if no other status can be identified
References the Medication resource.
References the Patient resource.
References the Observation resource. Information on the Most Recent Height and Weight observations will be shared.
Date the Medication was prescribed
References the Practitioner resource.
SIG
Medication supply information

Description
This will include the entire MEDENT medication SIG.

Description
Prescribed-Discontinued date period
Number of refills indicated for the medication order.
Quantity prescribed
Expected Discontinue Days
References the Organization resource.