

(g)(10) Standardized API - MEDENT VERSION 23.7**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

MedicationRequest

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the Medication
status	code	TRUE	Value will indicate if intended (Prescribed Date AFTER current date), active, or completed (Discontinued Date is populated and prior to current date or is current date)
medication	Reference	TRUE	References the Medication resource.
subject	Reference	TRUE	References the Patient resource.
effectivePeriod	Period	FALSE	The Prescribed Date (Start) and Discontinued Date (End) if populated
dosage	Dosage	FALSE	SIG

Dosage

Name	Type	Is Optional?	Description
text	string		This will include the entire MEDENT medication SIG.

This currently only includes items marked as Other Prescriber in the medication record.
While not part of the Core6.1.0 Requirements, this is shared to ensure data is available.