

(g)(10) Standardized API - MEDENT VERSION 23.7

Bundle

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Observation

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the Survey Record. This will reflect the Document UUID with ".hospice" appended
effective	dateTime	TRUE	Timestamp of the observation
status	code	FALSE	Value will be "final"
category	CodeableConcept	FALSE	Value will be "survey"
code	CodeableConcept	FALSE	Value will be LOINC 45755-6 Hospice Care [Minimum Data Set]
subject	Reference	FALSE	References the Patient resource.
valueCodeableConcept	CodeableConcept	FALSE	This will a SNOMED value for the checked checkbox value: Yes - 373066001 No - 373067005

At this time, this will be used to support additional documentation for: Hospice Care Screenings

Representation is based on the Hospice Care Guidelines from the PACIO Advance Healthcare Directive Interoperability IG

[ADI PMO Hospice Observation - PACIO Advance Healthcare Directive Interoperability Implementation Guide v2.0.0-ballot](#)