

(g)(10) Standardized API - MEDENT VERSION 23.5 R4**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Observation

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the lab test
status	code	FALSE	Value will be "final" if the lab has been marked Acquired or "registered" if there is no acquired date
category*	CodeableConcept	FALSE	Value will be "laboratory"
code	CodeableConcept	FALSE	LOINC code associated with the lab test
subject	Reference	FALSE	References the Patient resource.
performer*	Reference	TRUE	References the Practitioner resource.
effectiveDateTime	dateTime	TRUE	Date/time when the Lab Order was acquired.
valueQuantity	Quantity	TRUE	Lab result and indicated units. System will be from http://unitsofmeasure.org
note*	Annotation	TRUE	Note by Provider if available
hasMember*	Reference	TRUE	References the Observation information. If a lab test is part of a Panel, this will contain an array of references to each observation in the Panel. If the lab is not part of a Panel, this node will be omitted.