

(g)(10) Standardized API - MEDENT VERSION 23.5 R4

Bundle

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Observation

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the Social History record
effective	dateTime	TRUE	Timestamp of the observation
status	code	FALSE	Value will be "final"
category	CodeableConcept	TRUE	Value will be "social-history"
code	CodeableConcept	FALSE	Value will be the tobacco use status
subject	Reference	TRUE	References the Patient resource.
issued	instant	FALSE	Date of the observation
performer*	Reference	TRUE	References the Practitioner resource.
valueCodeableConcept	CodeableConcept	FALSE	Sent if the smoking status is based on an ICD10 code.
valueBoolean	boolean	TRUE	Sent if the smoking status is documented anywhere else. Will indicate true if the patient is a tobacco user.