

(g)(10) Standardized API - MEDENT VERSION 23.5 R4

Bundle

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Observation

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the vital entry
effective	dateTime	TRUE	Date/Time of the vital entry
status	code	FALSE	Value will be "final"
category*	CodeableConcept	TRUE	Value will be "vital-signs".
code	CodeableConcept	FALSE	LOINC of the vital entry
subject	Reference	FALSE	References the Patient resource.
performer*	Reference	TRUE	References the Practitioner resource.
valueComponent	Component	TRUE	Sent for Blood Pressure readings.
valueQuantity	Quantity	TRUE	Sent for all other vitals.
bodySite	CodeableConcept	TRUE	Body site associated with the vital