

**(g)(10) Standardized API - MEDENT VERSION 23.7****Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

**Observation**

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the lab test
status	code	FALSE	Value will be "final" if the lab has been marked Acquired or "registered" if there is no acquired date
category*	CodeableConcept	FALSE	Value will be "laboratory"
code	CodeableConcept	FALSE	LOINC code associated with the lab test
subject	Reference	FALSE	References the Patient resource.
valueQuantity	Quantity	TRUE	Lab result and indicated units. System will be from <a href="http://unitsofmeasure.org">http://unitsofmeasure.org</a>
hasMember*	Reference	TRUE	References the Observation information. If a lab test is part of a Panel, this will contain an array of references to each observation in the Panel. If the lab is not part of a Panel, this node will be omitted.
specimen	Reference	TRUE	References the Specimen resource.