

(g)(10) Standardized API - MEDENT VERSION 23.7

Bundle

Name	Type	Is Optional?
id	string	TRUE
meta	instant	TRUE
type	code	FALSE
total	int64	TRUE
link*	url	TRUE

Observation

Name	Type	Is Optional?
identifier*	Identifier	TRUE
effective	Period	TRUE
status	code	FALSE
category	CodeableConcept	TRUE
code	CodeableConcept	FALSE
subject	Reference	FALSE
valueCodeableConcept	CodeableConcept	FALSE
component:Industry	CodeableConcept	TRUE

Description
processId-GeneratedDateTime in UTC format
UTC format timestamp of the lastUpdate instant.
searchset
Total number of resources included in the response.
Link for the included data

Description
UUID of the Social History record
Start and/or Stop Date from Employment Information
Value will be "final"
Value will be "social-history"
History of occupation
References the Patient resource.
Occupation value from Occupation ONETSOC Detail (ODH)
Industry value based on Industry NAICS Detail (ODH) or free text value