

**(g)(10) Standardized API - MEDENT VERSION 23.7**

**Bundle**

Name	Type	Is Optional?
id	string	TRUE
meta	instant	TRUE
type	code	FALSE
total	int64	TRUE
link*	url	TRUE

**Observation**

Name	Type	Is Optional?
identifier*	Identifier	TRUE
effective	dateTime	TRUE
status	code	FALSE
category	CodeableConcept	TRUE
code	CodeableConcept	FALSE
subject	Reference	FALSE
valueCodeableConcept	CodeableConcept	FALSE
valueString	String	TRUE

Description
processId-GeneratedDateTime in UTC format
UTC format timestamp of the lastUpdate instant.
searchset
Total number of resources included in the response.
Link for the included data

Description
UUID of the Social History record
Timestamp of the observation
Value will be "final"
Value will be "social-history"
Value will be the LOINC for the Screening Question ex. if screening is "Do you want to talk about contraception or pregnancy prevention during your visit today" - LOINC 98076-3 will be sent
OR
Intent code if not tied to a Screening - ex. History Component selection
References the Patient resource.
Value will be sent based on the response to the screening. This will either contain SNOMED Intent Codes or the LOINC for the Screening Response ex. "code": "LA26438-4", "display": "Yes, I want to become pregnant", "system": "http://loinc.org"
Sent if the code is the Intent Code to display SNOMED or text description of Pregnancy Intent