

(g)(10) Standardized API - MEDENT VERSION 23.7

Bundle

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Observation

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the lab test
status	code	FALSE	Value will be "final" for shared screening documents and responses
category*	CodeableConcept	FALSE	value will be "survey"
code	CodeableConcept	FALSE	Question LOINC information
subject	Reference	FALSE	References the Patient resource.
performer*	Reference	TRUE	References the Practitioner resource.
effectiveDateTime	dateTime	TRUE	Document Date
valueCodeableConcept	CodeableConcept	TRUE	Answer LOINC Information
derivedFrom	Reference	TRUE	References the QuestionnaireResponse