

(g)(10) Standardized API - MEDENT VERSION 23.7

Bundle

Name	Type	Is Optional?
id	string	TRUE
meta	instant	TRUE
type	code	FALSE
total	int64	TRUE
link*	url	TRUE

Observation

Name	Type	Is Optional?
status	code	FALSE
category	CodeableConcept	TRUE
code	CodeableConcept	FALSE
subject	Reference	FALSE
valueCodeableConcept	CodeableConcept	TRUE

Description
processId-GeneratedDateTime in UTC format
UTC format timestamp of the lastUpdate instant.
searchset
Total number of resources included in the response.
Link for the included data

Description
Value will be "final"
Value will be "social-history"
Value will be the LOINC 76690-7
References the Patient resource.
SNOMED code of the Sexual Orientation value selected if any