

**(g)(10) Standardized API - MEDENT VERSION 23.7**

**Bundle**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>
id	string	TRUE
meta	instant	TRUE
type	code	FALSE
total	int64	TRUE
link*	url	TRUE

**Observation**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>
identifier*	Identifier	TRUE
effective	dateTime	TRUE
status	code	FALSE
category	CodeableConcept	TRUE
code	CodeableConcept	FALSE
subject	Reference	TRUE
issued	instant	FALSE
performer*	Reference	TRUE
valueCodeableConcept	CodeableConcept	FALSE
valueBoolean	boolean	TRUE

Description
processId-GeneratedDateTime in UTC format
UTC format timestamp of the lastUpdate instant.
searchset
Total number of resources included in the response.
Link for the included data

Description
UUID of the Social History record
Timestamp of the observation
Value will be "final"
Value will be "social-history"
Value will be the tobacco use status
References the Patient resource.
Date of the observation
References the Practitioner resource.
Sent if the smoking status is based on an ICD10 code.
Sent if the smoking status is documented anywhere else. Will indicate true if the patient is a tobacco user.