

(g)(10) Standardized API - MEDENT VERSION 23.7

Bundle

Name	Type	Is Optional?
id	string	TRUE
meta	instant	TRUE
type	code	FALSE
total	int64	TRUE
link*	url	TRUE

Observation

Name	Type	Is Optional?
identifier*	Identifier	TRUE
effective	dateTime	TRUE
status	code	FALSE
category*	CodeableConcept	TRUE
code	CodeableConcept	FALSE
subject	Reference	FALSE
performer*	Reference	TRUE
valueComponent	Component	TRUE
valueQuantity	Quantity	TRUE
bodySite	CodeableConcept	TRUE

Description
processId-GeneratedDateTime in UTC format
UTC format timestamp of the lastUpdate instant.
searchset
Total number of resources included in the response.
Link for the included data

Description
UUID of the vital entry
Date/Time of the vital entry
Value will be "final"
Value will be "vital-signs".
LOINC of the vital entry
References the Patient resource.
References the Practitioner resource.
Sent for Blood Pressure readings and Pulse Oximetry related readings for Inhaled O2 Concentration and Flow Rate
Sent for all other vitals.
Body site associated with the vital