

(g)(10) Standardized API - MEDENT VERSION 23.5 R4**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Organization

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	Practice ID
active	boolean	FALSE	Indicates if the Organization is active
type*	CodeableConcept	TRUE	Value will be "Healthcare Provider"
name	string	FALSE	The Patient Home, Work, or Mobile number or Email.
telecom*	ContactPoint	TRUE	Telecom information
address*	Address	TRUE	Address information
contact*	contact	TRUE	Practice contact information

Address

Name	Type	Is Optional?	Description
use	string	TRUE	Value will be "work"
line*	string	TRUE	Address Line
type	string	TRUE	Value will be "postal"
address-city	string	TRUE	City
address-state	string	TRUE	State
address-postalCode	string	TRUE	Postal Code
address-country	string	TRUE	Country
address-period	Period	TRUE	Period address is in use

ContactPoint

Name	Type	Is Optional?	Description
system	string	TRUE	Value will be phone or email.
value	string	TRUE	Phone or Fax #
use	string	TRUE	Value will be work.

Contact

Name	Type	Is Optional?	Description
name	HumanName	TRUE	Name of the Contact
telecom	ContactPoint	TRUE	Telecom information
address	Address	TRUE	Address information