

(g)(10) Standardized API - MEDENT VERSION 23.5 R4

Bundle

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Patient

Name	Type	Is Optional?	Description
identifier*	Identifier	FALSE	UUID of the Patient Record
active	Boolean	TRUE	Indicates if the Patient is active
name	HumanName	FALSE	Name
telecom*	ContactPoint	TRUE	The Patient Home, Work, or Mobile number or Email.
gender	string	FALSE	Administrative Gender of the Patient
birthdate	date	TRUE	Date of Birth
deceased	boolean dateTime	TRUE	deceased dateTime will be included if available. If not, the boolean will be included.
address*	Address	TRUE	Address information
maritalStatus	CodeableConcept	TRUE	Marital Status
contact*	Contact	TRUE	HIPAA Contact information
communication*	Communication	TRUE	Documented language preference
generalPractitioner*	Reference	TRUE	References the Practitioner resource

managingOrganization	Reference	TRUE	References the Organization resource
extension (race)	Extension	TRUE	Race
extension (ethnicity)	Extension	TRUE	Ethnicity
extension (birthsex)	Extension	TRUE	Sex at birth

Address

Name	Type	Is Optional?	Description
use	string	TRUE	Value will be "home"
line*	string	TRUE	Address Line
type	string	TRUE	Value will be "postal"
address-city	string	TRUE	City
address-state	string	TRUE	State
address-postalCode	string	TRUE	Postal Code
address-country	string	TRUE	Country
address-period	Period	TRUE	Period address is in use

HumanName

Name	Type	Is Optional?	Description
use	string	TRUE	Value will be "usual"
family	string	TRUE	Family name
given*	string	TRUE	Given name
suffix*	string	TRUE	Suffix

ContactPoint

Name	Type	Is Optional?	Description
system	string	TRUE	Value will be phone or email.
value	string	TRUE	Phone or email
use	string	TRUE	Value will be home or mobile

Contact

Name	Type	Is Optional?	Description
relationship*	CodeableConcept	TRUE	Relationship to the Patient
name	HumanName	TRUE	Name of the Patient Contact

telecom*	ContactPoint	TRUE	Telecom information
address*	Address	TRUE	Address information
period	Period	TRUE	Period address is in use

Communication

Name	Type	Is Optional?	Description
preferred	boolean	FALSE	Preferred language for communication
language	CodeableConcept	TRUE	Language indicated in preference