

(g)(10) Standardized API - MEDENT VERSION 23.7**Bundle**

Name	Type	Is Optional?
id	string	TRUE
meta	instant	TRUE
type	code	FALSE
total	int64	TRUE
link*	url	TRUE

Patient

Name	Type	Is Optional?
identifier*	Identifier	FALSE
active	Boolean	TRUE
name	HumanName	FALSE
telecom*	ContactPoint	TRUE
gender	string	FALSE
birthdate	date	TRUE
deceased	boolean dateTime	TRUE
address*	Address	TRUE
maritalStatus	CodeableConcept	TRUE
contact*	Contact	TRUE
communication*	Communication	TRUE
generalPractitioner*	Reference	TRUE
managingOrganization	Reference	TRUE
race	Extension	TRUE
ethnicity	Extension	TRUE
tribalAffiliation	Extension	TRUE
birthsex	Extension	TRUE

sex	Extension	TRUE
genderIdentity	Extension	TRUE
gender	Extension	TRUE

Address

Name	Type	Is Optional?
use	string	TRUE
line*	string	TRUE
type	string	TRUE
address-city	string	TRUE
address-state	string	TRUE
address-postalCode	string	TRUE
address-country	string	TRUE
address-period	Period	TRUE

HumanName

Name	Type	Is Optional?
use	string	TRUE
family	string	TRUE
given*	string	TRUE
suffix*	string	TRUE

ContactPoint

Name	Type	Is Optional?
system	string	TRUE
value	string	TRUE
use	string	TRUE

Contact

Name	Type	Is Optional?
relationship*	CodeableConcept	TRUE
name	HumanName	TRUE
telecom*	ContactPoint	TRUE
address*	Address	TRUE
period	Period	TRUE

Communication

Name	Type	Is Optional?
preferred	boolean	FALSE
language	CodeableConcept	TRUE

Description
processId-GeneratedDateTime in UTC format
UTC format timestamp of the lastUpdate instant.
searchset
Total number of resources included in the response.
Link for the included data

Description
This will include the following values and identifier-uses. Resource transactions should reference the official UUID identifier: usual: MRN official: UUID secondary: Account Number
Indicates if the Patient is active
Name
The Patient Home, Work, or Mobile number or Email.
Administrative Gender of the Patient
Date of Birth
deceased dateTime will be included if available. If not, the boolean will be included.
Address information
Marital Status
HIPAA Contact information
Documented language preference
References the Practitioner resource
References the Organization resource
Race
Ethnicity
Tribal Affiliation in Patient Information
Sex at birth

Sex field - supported values are male female other unknown asked-declined
GenderIdentity in Patient Information
Sex field - supported values are male female other unknown asked-declined

Description
Value will be "home"
Address Line
Value will be "postal"
City
State
Postal Code
Country
Period address is in use

Description
Value will be "usual"
Family name
Given name
Suffix

Description
Value will be phone or email.
Phone or email
Value will be home or mobile

Description
Relationship to the Patient
Name of the Patient Contact
Telecom information
Address information
Period address is in use

Description
Preferred language for communication
Language indicated in preference