

(g)(10) Standardized API - MEDENT VERSION 23.5 R4**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Practitioner

Name	Type	Is Optional?	Description
identifier	Identifier	FALSE	Identifier will be the provider NPI if available.
active	Boolean	TRUE	Indicates if Provider is active
name	HumanName	FALSE	Name of the Provider
telecom*	ContactPoint	TRUE	The Provider Home, Work, or Mobile number or Email.
address*	address	TRUE	Provider Address

HumanName

Name	Type	Is Optional?	Description
use	string	TRUE	Value will be "usual"
family	string	TRUE	Family name
given*	string	TRUE	Given name
suffix*	string	TRUE	Suffix

ContactPoint

Name	Type	Is Optional?	Description
system	string	TRUE	Values may include "phone" or "email"
value	string	TRUE	Phone or email
use	string	TRUE	Values may include "home," "work," or "mobile"

Address

Name	Type	Is Optional?	Description
use	string	TRUE	Value will be "work"
line*	string	TRUE	Address Line
type	string	TRUE	Value will be "postal"
address-city	string	TRUE	City
address-state	string	TRUE	State
address-postalCode	string	TRUE	Postal Code

address-country	string	TRUE	Country
address-period	Period	TRUE	Period address is in use